



Monitors & Signal Unit

SELECTION CHART

Use this feature comparison chart to select the device that's right for you

Fall Alert Wireless Signal Unit BPRC-64EA2
Bed Exit Monitor M-1000
Mobile Call Monitor MC-700
Mobile Monitor MM-600
Fall Monitor FM-400
Personal Monitor PM-800
Bed Tender® A14060900
Chair Tender® A14070900
Nurse Call Interface IS-200
Fall Alarm Interface IS-100
Chair Alarm Interface IS-600
Portable Nurse Call IS-300

	Fall Alert Wireless Signal Unit	Bed Exit Monitor	Mobile Call Monitor	Mobile Monitor	Fall Monitor	Personal Monitor	Bed Tender®	Chair Tender®	Nurse Call Interface	Fall Alarm Interface	Chair Alarm Interface	Portable Nurse Call
Dedicated System Specifically for alerting nurse's station of potential fall risk	●											
Automatic Reset Greatly reduces false alarms caused by normal movements and ensures device is always on	●	●	●	●	●				●	●	●	
AC Powered Provides reliable power without the use of batteries		●					●					
Battery Powered Allows for mobility; use single device for bed/chair settings; eliminates electrical cords	●		●	●	●	●		●	●	●	●	●
Incontinence Detection Immediately alerts staff of incontinence so skin deterioration can be prevented							●	●				
Nurse Call System Connection Alerts nurse's station of potential fall risk		●	●				●		●	●		
Fixed Installation Lets caregivers focus on providing care, instead of asset tracking									●	●		
Variable Delay Reduces false alarms by allowing for brief patient movement before sounding alert	●	●	●				●					
No Sensor Pad Required An economical monitoring solution effective in a variety of settings						●						●